

#### Command 0 - All Clear (Reset)

{0}

#### Command 1 - Channel Setup

{ 1, Channel, Operation, Post-Processing, FFT Samples }

\* Parameter value marked with asterisk are initial defaults.

Channel	Operation	Post-Processing	FFT Samples
0 Clear all channels	-	-	-
*1 Channel 1 2 Channel 2 3 Channel 3	0 Clear the selected channel. *1 Auto-ID 2 Voltage ( $\pm 10V$ ) (for Voltage probe) 4 Resistance 5 Period 6 Frequency 7 Temperature (Celsius) 8 Temperature (Fahrenheit) 9 Light 10 Voltage (0-5V) 11 Unavailable	*0 None 1 d/dt 2 d/dt, d <sup>2</sup> /dt <sup>2</sup> 10 FFT-Real 11 FFT-Real, FFT-Imaginary	- 1 to 13 : Samples used (*6) 2^n (2-8192)
4			
5			
6 Unavailable			
10			
11			
12			
13 Built-in accelerometer channel (X axis)			
14 Built-in accelerometer channel (Y axis)			
15 Built-in accelerometer channel (Z axis)			

#### Command 3 - Sample and Trigger Setup

{ 3, Sampling Interval, Number of Samples, Record Time, Trigger Source, Trigger Threshold, Trigger Edge, Clock Source }

Sampling Interval	Number of Samples	Record Time	Trigger Source	Trigger Threshold	Trigger Edge	Clock Source
0.00002	Number of seconds	1-120000	Number of samples	0 Off	-	-
- 16000		(*100)		*1 Absolute time recording		
(*0.1)				2 Relative time recording		
				None 0 Start sampling immediately upon receipt of Command 3.	-	-
				Key Trigger *1 Start measuring when power on/off button is pressed once.		
				Threshold	0 Falling edge	*0 Timer (Sampling interval)
				2 CH1 Trigger 3 CH2 Trigger 4 CH3 Trigger 5 CH1 Trigger 6 CH2 Trigger 7 CH3 Trigger	*1 Rising edge	10 Same as Trigger Source (1-5)
				8 9 10 Unavailable 11	2 Rising and falling edge	
				13 Built-in accelerometer channel (X axis) 14 Built-in accelerometer channel (Y axis) 15 Built-in accelerometer channel (Z axis)	Count Number (sec)	
				20 Countdown Trigger	1 - 10 (*10)	
				-1 Command 8		

#### Command 4 - Conversion Equation Setup

{ 4, Equation Number, Equation Type, Number Format, Constants }

Equation Number	Equation Type	Number Format	Constants
*0 Clear all equations.	-	-	-
1 Equation 1 (Channel 1)	*0 Clear equation selected by the equation number parameter.	*0 Standard	K0(, K1, ..., K9) *3,4
2 Equation 2 (Channel 2)	1 Polynomial K0+K1X+K2X <sup>2</sup> +...+KnXn *1	10 Integer part (Decimal part cut off)	K-4(, ..., K-1, K0, K1, ..., K5) *3,4
3 Equation 3 (Channel 3)	2 Mixed polynomial K-mX-m+...+K-1X-1+K0+K1X+...+KnXn *2		
	3 Power K0K(X)+K2		
	4 Modified power K0K1(X)+K2		
	5 Logarithmic K0+K1 ln(X)		
	6 Modified logarithmic K0+K1 ln(1/X)		
	7 Exponential K0 e(K1X)+K2		
	8 Modified exponential K0 X(K1X)+K2		
	9 Geometric K0 X(K1X)+K2		
	10 Modified geometric [K0+K1 ln(K2X)]-1+K3		
	11 Reciprocal logarithmic [K0+K1 ln(K2X)]-1+K3		
	12 Steinhart-Hart mode [K0+K1 ln(1000X)+K2(ln 1000X)3]-1+ K3		
13 Equation 4 (Built-in accelerometer channel (X axis))			
14 Equation 5 (Built-in accelerometer channel (Y axis))			
15 Equation 6 (Built-in accelerometer channel (Z axis))			
4 Unavailable			

\*1 Polynomial: Input constants in sequence, from n = 0 to 9.

\*2 Mixed polynomial: Input constants in sequence from m = 4 to 1, and n = 0 to 5.

\*3 Input of zero for constants can be skipped if all remaining constants are not used.

\*4 Input 0 for constants that are not used.

\* When the conversion result of the "conversion equation" selected by Command 4 causes an overflow, data-logger sends a result of zero (0) to the calculator.

#### Command 5 - Data Range Setup

{ 5, Channel Select, Data Select, Data Begin, Data End, Step, K, FFT Samples }

Channel Select	Data select	Data Begin	Data End	Step	K	FFT Samples
*0 Current send channel	*0 Raw data					
1 Channel 1	1 d/dt					
2 Channel 2	2 d <sup>2</sup> /dt <sup>2</sup>					
3 Channel 3	9 A/D conversion value	1-120000	1-120000	Data Range Steps (*1)		
	10 FFT-Real	(*1)	(*1)	-1: Data range number / K (*1)	+255	1-13 : Samples used 2^n (2-8192) (*6)
	11 FFT-Imaginary					
4 Unavailable						
5 Recorded time data						
10 Unavailable						
13 Built-in accelerometer channel (X axis)						
14 Built-in accelerometer channel (Y axis)						
15 Built-in accelerometer channel (Z axis)						

#### Command 6 - System Setup

{ 6, Command, Auto Power Off Time }

Command
0 Abort Sampling
1 Unavailable
2
3 Turns sound off
4 Turns sound on
10 Unavailable

#### Command 7 - Status Check

{ 7 }

#### Command 8 - Sampling Start

{ 8 }

#### Command 10 - Sensor Warm-up

{ 10, Warm-up Time (sec) }

Warm-up Time (sec)
0.1 - 360
0 Auto
-1 None
-2 Normal warm-up

#### Command 11 - Buzzer and LED Operation Commands

{ 11, Output Select, Length, Period }

Output Select	Length	Period
*0 Buzzer		Sound period (100μsec unit)
1 Unavailable		
2 Ready LED		
3 Sampling LED		
4 Error LED		
5 Battery LED		
	Operation Time (100μsec unit)	

#### Command 12 - Data Send Sequence

{ 12, Send Sequence }

Send Sequence
*0 Non-real Time Format
1 Real Time Format

#### Command 21 - Auto ID Information detect

{ 21 }